

The Bullard™ Laryngoscopes

A faster, less traumatic solution
for the difficult intubation.



Designed for difficult
airways yet versatile enough
for routine intubation.



For Difficult and Routine Intubations

The Adult, Child and Newborn/Infant Bullard Elite™ Laryngoscopes are anatomically designed to facilitate rapid visualization of the larynx without requiring alignment of the pharyngeal, laryngeal and oral axes.

The Bullard Laryngoscope was developed for the difficult-to-intubate patient but is simple enough for routine intubation of patients with normal airways. Members of the anesthetic care team typically attain proficiency within two to four cases.

Bullard Elite Laryngoscopes can assist you in managing difficult, life-threatening conditions such as:

- Cervical spine fracture
- Upper body burns
- Trauma
- Patients with TMJ immobility
- Congenital airway problems such as Pierre Robin syndrome



Laryngoscope Kits*

Catalog No.	Description
LARS-AL	Adult Bullard Elite™ Laryngoscope w/Multifunctional Stylet
LARS-A	Adult Bullard Elite™ Laryngoscope w/Standard Introducing Stylet
LARS-PL	Child (Pediatric Long) Bullard Elite™ Laryngoscope w/Std. Introducing Stylet
LARS-P	Newborn/Infant (Pediatric) Bullard Elite™ Laryngoscope, no stylet included

Intubating Stylets

Catalog No.	Description
LMS-AL	Multifunctional Stylet, Adult Large, 5.6mm OD / 4.5mm ID
LIS-A	Standard Introducing Stylet, Adult
LIS-PL	Standard Introducing Stylet, Child (Pediatric Long)

Optional Accessories

Catalog No.	Description
LAR-AE†	Blade Extender for Adult Bullard™ Laryngoscope (6/pkg)
LAR-ER	Bullard™ Laryngoscope Blade Extender Remover
LBH-1	Bullard™ Laryngoscope Battery Handle
LFK-3	Bullard Elite™ Laryngoscope Fiberoptic Adapter
ALU-2B	Fiberoptic Light Source, dual lamp (150 W, 115 V, 60 Hz)
G93	High Transmission Autoclavable Fiberoptic Light Guide, "snap-on," 4mm, 7ft.
98-0340	Adapter to Storz Light Guides
98-0341	Adapter to Wolf Light Guides
8914	Adapter to connect ACMI Light Guides to Olympus Light Sources
75U	Universal Stopcock, LeurLok®
AE-1	Endoscopic Instrument Lubricant (3/pkg)

Replacement Parts

Catalog No.	Description
LB-300	Cleaning Brush
LSL-3E	Spare Lamp Package (3/pkg)
LAS-B	Anti-Fog Solution

† Supplied sterile, single-use disposable.

*All Bullard Elite Laryngoscope Kits include the following accessories: Cleaning Brush, Anti-Fog Solution, Spare Lamp and Carrying Case.

For more information,
please contact your local GYRUS ACMI
sales representative or call
Customer Service at 888.524.7266


GYRUS ACMI

136 Turnpike Road
Southborough, MA 01772 USA

888.524.7266 PHONE
262.639.7205
262.681.4038 FAX
www.acmicorp.com WEB

© 2005 ACMI Corporation. All rights reserved. ™ and ® Trademarks and Registered Trademark, respectively, of ACMI Corporation or its Subsidiaries in the U.S. and/or other countries of the world. For U.S. Distribution Only.

SS501-1105 Subject to change without notice